

# C.BIRD™ | An Innovative Cell Culture for Hybridoma in Early Stages

Yun Chen, MS, Charles Tsai, PhD  
CYTENA Bioprocess Solutions, Taipei City, Taiwan



## Abstract

Hybridoma technology is the first and well-established methodology to produce monoclonal antibodies (mAbs). However, it faces the obstacles of the monoclonal cell line development, such as low growth rate and poor productivity under the low initial density. By growing single colonies in microwells, the cells endure an extended time in static culture, in which the concentration of nutrients and oxygen in the medium decreases over time. To optimize the clonal expansion from low cell density and to speed-up the antibody production for the early-stage analysis, we

introduce the C.BIRD™, our innovative technology. C.BIRD™ offers continuous mixing in 96-well plates to improve the culture environment for hybridoma cells by providing a higher diffusion rate of oxygen for a more homogeneous culture environment. In this application note, we compared the cell growth and mAb production efficiencies under different culture conditions. As a result, we observed increased cell growth and significantly improved mAb production in C.BIRD™ cultures in comparison to static or mixing cultures. Our C.BIRD™ technology holds a great potential to improve the hybridoma development by shortening the timeline of culturing and early-stage analytical studies.

## Introduction

Monoclonal antibodies (mAbs) have been playing a fundamental role in biomedical industries since their discovery. The high specificity, affinity, and homogeneous nature of mAbs make them valuable in various research techniques, including flow cytometry, immunoassays, diagnostics, and also as therapeutics. During the development of mAbs, the antigen of interest is injected in a mouse and the spleen cells are isolated. After fusing the spleen cells and myeloma cells, hybridomas were distributed into microwell plates with low cell density to develop single colonies. For the humanized mAb development, mAbs should be altered to meet the requirement of human immune system in which genetically-engineered hybridomas go through another clonal selection. The time-consuming process of selecting and expanding desirable clones is one of the obstacles, slowing down the development of hybridoma workflow.

Due to the lack of nutrients and oxygen supplements during the long culture period, the poor growth rate of the cells is one of the major causes for the longdeveloping process of the technology.

In this study, we proposed a new method to distribute the nutrients and increase the oxygen diffusion in the culture medium by using the C.BIRD™. The continuous mixing of the C.BIRD™ provides a homogeneous suspension culture in microwell plates. In our result, hybridomas possessed a higher growth rate and mAb productivity by using the C.BIRD™ compared to the conventional culture methods. Undoubtedly, our method shortens the cell culture expansion process, which starts from a single colonie or a low cell density cell culture. Therefore, the C.BIRD™ has the potential to facilitate the early-stage analytic screening by shortening the required timeline for early mAb development.




Figure 1. Diagram of experiment design.

## Materials and methods

A mouse mAb-expressing hybridoma cell line (inhouse) was used in this study. Cells were cultured in Dulbecco's Modified Eagle Medium, High glucose (#CC103-0500, GenedireX) supplemented with 10% FBS (Gibco) and 1% penicillin/streptomycin. Standard 96-well plates (#0030730011, Eppendorf) were used for C.BIRD™ and static cultures (Figure 1). Flasks (#431143, Corning) were placed on a shaker for mixing cultures (shaking speed: 90 rpm; orbital: 20mm). Cultures were kept at 37°C in a water-jacketed incubator with 5% CO<sub>2</sub>.


The initial cell concentration was  $2.5 \times 10^3$  cells/ml in all cell culture conditions. The volume of the medium in both static and the C.BIRD™ groups was 200 µl/well in standard 96-well plates, while the volume of the shaker-flask culture was 30 ml in total. Total and viable cells number and viability were counted with an automatic cell counter (TC20 from Bio-Rad) from Day 4 to Day 7. To compare the antibody production efficiencies in each cell culture method, we measured the antibody concentration using mouse IgG (Total) Uncoated ELISA kit (# 88-50400-88, Invitrogen) in the supernatants of all cultures collected at day 7. Data were analyzed by one-way analysis of variance (ANOVA). Significance of p-value is listed as the following: 0.12(ns), 0.033(\*), 0.002(\*\*), 0.0002(\*\*\*) and <0.0001(\*\*\*\*). Data are shown as mean ± SD.

## Results and discussion

The C.BIRD™ provides a perpendicular mixing in the multi-well plates which introduces a better circulation of the nutrients and oxygen by the constant flow. To compare whether the mixing of the C.BIRD™ method could improve the cell culture conditions of hybridoma, we designed the experiment as shown in Figure 1. We measured total cell density, viable cell density, and analyzed the viability of the three different cell culture methods. Under continuous mixing (25s/cycle) in the C.BIRD™ culture in the standard 96-well plate, the total cell density peaked at  $2.44 \times 10^6$  cells/ml on average at Day 6, which were 1.99-fold and 1.89-fold higher compared to the static and shaker-flask groups, respectively (Figure 2A).

Additionally, viable cell density of the C.BIRD™ group an average of  $2.04 \times 10^6$  cells/ml at Day 6 and the other two groups had only an average of  $9.22 \times 10^5$  cells/ml and  $1.11 \times 10^6$  cells/ml in static and shaker-flask groups, respectively (Figure 2B). The results showed significant improvement in cell growth in the C.BIRD™ cultures. C.BIRD™ provided more efficient conditions to expand the hybridoma cells from a relatively low cell density.

On the other hand, the cell viability in the C.BIRD™ group was not compromised, regardless of the high cell density (Figure 2C). On Day 6, while sustaining the highest viable cell density within three groups, the C.BIRD™ group had 83% of viability on average, close to that of the shaker-flask group (87%) and was higher than the static group (75%). The result demonstrated that the C.BIRD™ method could increase environmental carrying capability and helped to expand the cell population. To investigate the efficiency of Mouse mAb production, supernatants were collected at Day 7. We measured a high concentration of mAb in the C.BIRD™ group, which reached an average of 3309.2 µg/mL, and outcompeted that in both static and shaker-flask groups (Figure 2D).



**Figure 2.** The cell growth and antibody production of three culture methods in 7 days of continuous culture. The C.BIRD™ facilitated the cell growth in both **A)** total cell density and **B)** viable cell density while maintaining the same **C)** viability as the other two groups. In addition, the C.BIRD™ group showed a significantly higher **D)** antibody production at Day 7 compared to both static and shaker-flask groups.

## Conclusion

Hybridoma cultures share similar issues with low productivity and cell density, the C.BIRD™ method provided an innovative way to boost cell growth in an extremely low seeding cell density. Our study showed a magnificent improvement in 1) total cell density, 2) viable cell density, 3) maintenance of viability during long-term culture, and 4) high productivity of mAb, comparing to static culture in 200 µL/well conditions. More excitingly, the C.BIRD™ group outperformed the shaker-flask group in cell growth and antibody production. By using the C.BIRD™, researchers can accelerate clonal selection and shorten the timeline for hybridoma development.

## Reference

1. Lianchun F, Giovanni R, Kathleen B, Jun T, Joon Chong Y, Reb R, Tapan K D. Comparative study of therapeutic antibody candidates derived from mini-pool and clonal cell lines. *Biotechnol Prog.* 2017;33(6):1456-1462.  
[DOI: 10.1002/btpr.2477](https://doi.org/10.1002/btpr.2477).
2. Hilal Ahmed P, Shivangi S, Sweetie S, Tripti S, Shubbir A, Chandresh S, Rajesh K. Hybridoma technology a versatile method for isolation of monoclonal antibodies, its applicability across species, limitations, advancement and future perspectives. *Int Immunopharmacol.* 2020;85:106639.  
[DOI:10.1016/j.intimp.2020.106639](https://doi.org/10.1016/j.intimp.2020.106639).



©2023 BICO AB. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of BICO is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO may refer to the products or services offered by other companies by their brand name or company name solely for clarity and does not claim any rights to those third-party marks or names. BICO products may be covered by one or more patents. The use of products described herein is subject to BICO's terms and conditions of sale and such other terms that have been agreed to in writing between BICO and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of BICO products in practicing the methods set forth herein has not been validated by BICO, and such nonvalidated use is NOT COVERED BY BICO'S STANDARD WARRANTY, AND BICO HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner BICO's terms and conditions of sale for the instruments, consumables or software mentioned, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by BICO that it currently or will at any time in the future offer or in any way support any application set forth herein.

Edited version: August 2023 | CBS\_PUB\_CBIRD\_App-Note-05\_Digital