

# EZHold Lid Holder

The EZHold Lid Holders are accessories for the C.BIRD and S.NEST microbioreactors that keep lids from coming into direct contact with microwell plates. With the EZHold in place, you can easily add supplements or other components to wells while minimizing sample cross contamination.



## Features

- Separates lid and microwell plate
- Keeps fluidic tubes from coming into contact
- Reduces the risk of cross contamination between wells
- Easy-to-use design
- Compatible with 96- and 24-well formats
- Available in C.BIRD and S.NEST formats

## Specification

Dimension (WxDxH):  
17.8 x 22.3 x 20 cm

Weight:  
2.2 kg

## EZHold ordering information

EZHold ordering information		
Compatibility	C.BIRD	S.NEST
Item No.	AC0104	AS0101

## CONTACT

E-mail: [info@cytena-bps.com](mailto:info@cytena-bps.com)  
Website: [www.cytenua-bps.com](http://www.cytenua-bps.com)  
Taiwan phone: (+886) 2-2720-6135  
E.U. phone: (+46) 31-128-700  
U.S. phone: (+1) 833-235-5465

**CYTENA BPS**   
A BICO COMPANY

For research use only. Not for use in diagnostics.  
All product names and trademarks are the property of their  
respective owners. ©2021 BICO. All rights reserved.  
CBSFL001